

Samples for mineral oil analysis

Advice for sample shipment

How to handle the samples?

What kind of sample	Amount/size needed	Packaging for shipment
Food e.g. solid, liquid, powder   	min. 50g	glass jar/ PET cup with aluminium foil <div> <div> plastic lid aluminium foil glass jar </div>  <div> aluminium lid aluminium foil PET cup </div>  <div>Source: Packart.com GmbH</div> </div> <p>Please note: Make sure that the lids are tightly closed and avoid using adhesive tape for closing.</p> aluminium foil or polyamide bag <div>  <div> Polyamid bag (Rilsan® bag, material: PA 11) </div>  <div>Source: Tub-ex</div> </div>
Packaged Food 	min. 50g	aluminium foil 
Packaging material e.g. paper, cardboard, plastic foil, jute bag, big bag 	min. 50g minimum size for migration test: 8.3 inch x 11.7 inch 210 mm x 297 mm (5 sheets)	aluminium foil <div>     </div>

Which materials are not suitable as packaging materials?



plastic materials



jute



cardboard/paper